

Assemble Your Apollo Cart



STEP 1

Select a type 'E' or 'U'

Type E
60cm wide x 40cm deep



or


Type U
40cm wide x 60cm deep




STEP 2

Select a height for type E


4.6 Panels Cart
Height 209cm
APL-4.6E-TD/GD/SD/ND



4 Panels Cart
Height 184cm
APL-4.0E-TD/GD/SD/ND



3 Panels Cart
Height 146cm
APL-3.0E-TD/GD/SD/ND




2 Panels Cart
Height 109cm
APL-2.0E-TD/GD/SD/ND




STEP 2

Select a height for type U


4.6 Panels Cart
Height 209cm
APL-4.6U-TD/GD/SD/ND




4 Panels Cart
Height 184cm
APL-4.0U-TD/GD/SD/ND



3 Panels Cart
Height 146cm
APL-3.0U-TD/GD/SD/ND



2 Panels Cart
Height 109cm
APL-2.0U-TD/GD/SD/ND



STEP 3


Select wheels / levelling feet

	Description	Apollo Heights
	12.7cm wheel	reduce 0cm (default)
	10cm wheel	reduce 2cm
	level caster	reduce 7.5cm
	level feet	reduce 13cm


STEP 4

Select a cart width for E type


Single 75cm



Double 146cm




Triple 217cm




STEP 4

Select a cart width for U type


Single 55cm



Double 106cm




Triple 157cm




STEP 5

Choose drawer slider type


1. Tilt Panel
(up to 10kg per tray)



2. Combination of
Panels & Telescopic
Rails



3. Telescopic Rails
(up to 20kg per tray)



STEP 6

Choose a door type



solid door



glass door



tambour door



open cart
(no door)

STEP 7

Choose a locking system

1. No lock



2. Key lock



3. Electronic
code lock



4. RFID code
lock



STEP 8

Choose with drawer fronts or not



Choose a drawer configuration:

Configuration C1



Configuration C2



Configuration C3



STEP 9

Choose accessories

Aluminium Accessory Rail
E Type: #APL-SRAIL-E-KIT
U Type: #APL-SRAIL-U-KIT
For mounting various
accessories on carts.



Push Handle
#CRTA-ALPH-MK/NM
Aluminium 40cm



Disinfection Liquid Bottle Holder
#DLBH-PP +
#DIN-HKS-GBH-SET
Simple, hidden, click to
open mechanism
Max: 7cm bottle diameter



Trio Glove Holder
#TR-GBH-ABS +
#DIN-HKS-GBH-SET
Max box size 14.5 x 24.5 x 8cm
Requires 2 accessory rails



Glove Holder
#S-GBH-ABS +
#DIN-HKS-GBH-SET
Max box size 14.5 x 24.5 x 8cm



Waste Basket
#CRTA-WSBIN
Waste basket large – 10l with
Holder 26 x 20 x 36.5



STEP 9

Choose accessories

Foldable Frame
#APL-SFRAME-FL
For side shelf
Side tray



Side Basket
#CRTA-SBASKET
Wire side basket
44 x 16 x 19cm



Tilt Bin Holder 2
#APL-TBH-2
Tilt bin holder for 2 rows
Height above cart: 62.1 cm




Tilt Bin 6
#tiltbin6



Tilt Bin 5
#tiltbin5



Tilt Bin 4
#tiltbin4



STEP 9

Choose accessories for medication



Medication Cassette 3 Levels
ABS, cream
12 large drawers or 18 small
drawers

W595 x H339mm

MED-CAST-3L



Medication Drawer
ABS, cream
8 x 33.8cm – 80cm wide
MED-DWR-80



Medication Drawer
ABS, cream
12 x 33.8cm – 120cm wide
MED-DWR-120



Keylock provided
as standard



RFID electronic
lock with four
digit user code





info@distinctivemedical.com – 01928 571 801 – www.distinctivemedical.com