

INSAUSTI MEDICAL TROLLEY GENERAL CATALOGUE



Northern Ireland Framework Agreement – Clinical Trolleys and Equipment 6001970

Lot 2 – Resuscitation Trolleys

Distinctive Medical Products Ltd is pleased to have been awarded a position on the Northern Ireland Clinical Trolleys and Equipment Framework as an approved supplier for the provision of Insausti® resuscitation trolleys, associated accessories, and the wider Insausti® product range.

The trolley specifications detailed within this brochure represent a selection of our available solutions. Should these specifications not fully meet your clinical or operational requirements, our team would be delighted to discuss the development of a bespoke trolley configuration tailored to your specific needs.

Our experienced customer service team is available to provide guidance and support regarding trolley specifications, configurations, and any related product enquiries.

How to Place an Order

To place an order, please contact the Distinctive Medical Customer Service team by telephone or email, quoting the Northern Ireland Framework reference, product code, and required quantity.

Telephone: 01928 571801

Email: info@distinctivemedical.com

Please note that the standard lead time for all medical trolleys is approximately 28 days.


INSAUSTI TROLLEY TECHNICAL FEATURES

SQUARE TROLLEY

ISO 480 x 480 mm



SQS




SQM



DOUBLE TROLLEY

ISO 900 x 630 mm



300 SERIES




TECHNICAL FEATURES

SIMPLE TROLLEY

ISO 640 x 480 mm

ISO International Organization for Standardisation, ISO norm 3394
600 x 400 mm ISO drawers. Standard, interchangeable with hospital ISO furniture



LOW



300 SERIES



100 SERIES



TECHNICAL FEATURES

SQUARE TROLLEY

ISO
480 x 480 mm

SOFT CLOSING SLIDES (Optional)

No closing noise is produced Operating room, reanimation area, ICU



ASF



Prepared for elevated accessories

PE
SANITIZED



Ergonomic handle

Total opening
Steel ball-bearing slides with
retention system, 20 kg

PP
SANITIZED



Front stop
It protects the slide against
strong impacts when closing



TECHNICAL FEATURES

SIMPLE TROLLEY / DOUBLE TROLLEY

ISO

640 x 480 mm

900 x 630 mm

SOFT CLOSING SLIDES (Optional)

No closing noise is produced Operating room, reanimation area, ICU



Prepared for elevated accessories



ASF Simple trolley
ASFN Double trolley

PE
SANITIZED



Ergonomic handle

Total opening
Steel ball-bearing slides with
retention system, 20 kg

PE
SANITIZED



ABS

Front stop
It protects the slide against
strong impacts when closing

TECHNICAL FEATURES

MONOBLOC STRUCTURE

AISI-304 stainless steel

Assembly with rivets and TIG welding

Side panels prepared for accessories (holes with caps)



PE, PP SANITIZED

Withstand disinfectants for the elimination of pathogens

Withstand 100% pure alcohol (not diluted), bleach, detergent, soap solutions, oils, polish, aqueous ammonia, chlorhexidine

Withstand constant cleaning protocols, they do not degrade

PE, POLYETHYLENE WORK SURFACE

Protection on all 4 sides



TECHNICAL FEATURES

PP, POLYPROPYLENE SQUARE DRAWERS

ISO 480 x 480 mm

Polypropylene drawers, PP, manufactured in a single piece with slots for optional dividers.
Heights 57, 130 mm.



ISO DRAWERS

ABS ISO drawers of standard dimensions according to ISO 3394.
Heights 50, 100, 200 mm.

ISO 640 x 480 mm 900 x 630 mm



TECHNICAL FEATURES



STATIC TROLLEYS

SQS SERIES - 480 x 480 mm

- The work surface design allows you to move the trolley easily
- Ø 100 mm wheels, 2 antistatic, 2 with brake



300 SERIES - ISO 640 x 480 mm

- The work surface design allows you to move the trolley easily
- Ø 100 mm wheels, 2 antistatic, 2 with brake



300 SERIES - ISO 900 x 630 mm

- The work surface design allows you to move the trolley easily
- Ø 100 mm wheels, 2 antistatic, 2 with brake



TECHNICAL FEATURES

EXTENSIVE MOBILITY TROLLEYS

SQM SERIES - 480 x 480 mm

- Elevated push handle, left side
- Bumpers for hits
- Ø 125 mm simple wheels, 2 antistatic, 2 with brake



CE



I00 SERIES - ISO 640 x 480 mm

- Elevated push handle, left side
- Bumpers for hits
- Ø 125 mm double wheels, 2 antistatic, 2 with brake



CE



500 SERIES

- 480 x 480 mm
- 640 x 480 mm
- 900 x 570 mm

- Elevated push handle, left side
- Bumpers for hits
- Ø 125 mm simple wheels, 2 antistatic, 2 with brake



TECHNICAL FEATURES

SQUARE TROLLEY

Drawer colour



W White



B Blue



R Red

SIMPLE TROLLEY / DOUBLE TROLLEY

Front drawer colour



W White



B Blue



F Lilac



G Green



Y Yellow




R Red

TECHNICAL FEATURES

Plastic material. ISO norm 7810 identification card, international standard

SQUARE TROLLEY

Each drawer includes as standard one white card




White card

SIMPLE TROLLEY / DOUBLE TROLLEY

Each drawer includes as standard one card holder and one white card



Card holder



White card

COLOUR CARDS (Optional)



CR Red



CG Green



CP Pink



CB Blue



CO Orange



CF Lilac



CY Yellow



CN Grey



FRAMEWORK PRODUCT INFORMATION

LOT 2 - CORE LIST

EXTENSIVE MOBILITY RESUSCITATION TROLLEY

PPC: TDT002904

Product description:

- Dimensions
 - Width 640mm
 - Depth 480mm
 - Height 1015mm
- Monobloc structure of AISI-304 stainless steel
- Polyethylene work surface
- ISO drawers 600 x 400mm
 - 3 shallow
 - 2 deep
- Tray divider sets with four compartments
- Tray label holders and labels
- 125mm wheels (x2 antistatic x2 with brake)
- Tamper evident security seal
 - Black polyethylene bridles
 - Stainless steel bar
 - Disposable seal cutting system

Accessories:

- Gas bottle stand – left side
 - For 110/140mm bottles
 - Maximum 18kg
 - Adjustable black nylon strap with quick buckle locking
 - Stable bottle fixation for base of bottle
- Fixed IV pole – attached to the back of the trolley
 - Stainless steel
 - Height adjustable
 - 4kg load capacity
- ISO drawers 600 x 400mm
 - Available in ABS and/or transparent polycarbonate material
 - Sizes 8, 19 and 33 litres
 - All trays are removable for cleaning purposes
 - Tray runners preventing accidental withdrawal



LOW TROLLEYS, REDUCED HEIGHT

ISO

480 x 480 mm



STATIC TROLLEYS

To store it under 900 mm work tables

For users working in sitting position



Product code: Q001LW B, SQ061M4LW B.32R

Product description:

- DIMENSIONS 480 x 480 mm
- HEIGHT 855 mm (WORK SURFACE)
- WORK SURFACE, POLYETHYLENE PE
- MONOBLOC STRUCTURE, STAINLESS STEEL AISI-304
- SQUARE DRAWERS
- ø 100 mm DOUBLE WHEELS, 2 ANTISTATIC, 2 WITH BRAKE

Accessories

- 32R: SIDE SUPPORT



Product code: 3100LW W, 362M5LW W.32R

Product description:

- DIMENSIONS 640 x 480 mm
- HEIGHT 855 mm (WORK SURFACE)
- WORK SURFACE, POLYETHYLENE PE
- MONOBLOC STRUCTURE, STAINLESS STEEL AISI-304
- SQUARE DRAWERS
- ø 100 mm DOUBLE WHEELS, 2 ANTISTATIC, 2 WITH BRAKE

Accessories

- 32R: SIDE SUPPORT

MULTIFUNCTION SQS SERIES



STATIC TROLLEYS



Product code: Q025 W, SQ022 W.32R

Product description:

- DIMENSIONS 480 x 480 mm
- HEIGHT 980 mm (WORK SURFACE)
- WORK SURFACE, POLYETHYLENE PE
- MONOBLOC STRUCTURE, STAINLESS STEEL AISI-304
- SQUARE DRAWERS
- ø 100 mm DOUBLE WHEELS, 2 ANTISTATIC, 2 WITH BRAKE

Accessories

- 32R: SIDE SUPPORT

MULTIFUNCTION SQS SERIES

ISO

480 x 480 mm

ICU, CARDIOLOGY, DRESSING, OUTPATIENTS CLINICS



STATIC TROLLEYS

DIVIDERS (Optional)



8 DMQ, 8 DSQ



Product code: Q09I B, SQ06I B.26(RA4Q.A5).32R

Product description:

- DIMENSIONS 480 x 480 mm
- HEIGHT 980 mm (WORK SURFACE)
- WORK SURFACE, POLYETHYLENE PE
- MONOBLOC STRUCTURE, STAINLESS STEEL AISI-304
- SQUARE DRAWERS
- ø 100 mm DOUBLE WHEELS, 2 ANTISTATIC, 2 WITH BRAKE

Accessories

- 26(RA4Q.A5): BINS DISPENSERS
- 32R: SIDE SUPPORT

MULTIFUNCTION SQS SERIES

ISO
480 x 480 mm

ICU, CARDIOLOGY, DRESSING, OUTPATIENTS CLINICS



STATIC TROLLEYS

DIVIDERS (Optional)



8 DMQ, 8 DSQ



Product code: Q096 W, SQ06I W.26(RA4Q.A5).30.32L

Product description:

- DIMENSIONS 480 x 480 mm
- HEIGHT 980 mm (WORK SURFACE)
- WORK SURFACE, POLYETHYLENE PE
- MONOBLOC STRUCTURE, STAINLESS STEEL AISI-304
- SQUARE DRAWERS
- ø 100 mm DOUBLE WHEELS, 2 ANTISTATIC, 2 WITH BRAKE

Accessories

- 26(RA4Q.A5): BINS DISPENSERS
- 30: FOLD DOWN TRAY
- 32L: SIDE SUPPORT

MULTIFUNCTION SQS SERIES

ISO
480 x 480 mm

ICU, CARDIOLOGY, DRESSING, OUTPATIENTS CLINICS



STATIC TROLLEYS



DIVIDERS (Optional)



8 DMQ, 8 DSQ

Product code: Q092 W.(6)CP, SQ06I W.26(RA5Q.A5.RIQ).30.33L.(6)CP

Product description:

- DIMENSIONS 480 x 480 mm
- HEIGHT 980 mm (WORK SURFACE)
- WORK SURFACE, POLYETHYLENE PE
- MONOBLOC STRUCTURE, STAINLESS STEEL AISI-304
- SQUARE DRAWERS
- ø 100 mm DOUBLE WHEELS, 2 ANTISTATIC, 2 WITH BRAKE

Accessories:

- 26(RA5Q.A5.RIQ): BINS DISPENSERS – CENTRAL SUPPORT
- 30: FOLD DOWN TRAY
- 33L: DISPOSABLE CONTAINER SUPPORT
- (6)CP: PINK CARDS

MULTIFUNCTION 300 SERIES

ICU, CARDIOLOGY, DRESSING, OUTPATIENTS CLINICS

ISO
640 x 480 mm



STATIC TROLLEYS

DIVIDERS (Optional)



H5/8, (3)H10/4, H20/4



Product code: 315I B, 365 B.26LW(RI).33R.50

Product description:

- DIMENSIONS 640 x 480 mm
- HEIGHT 980 mm (WORK SURFACE)
- WORK SURFACE, POLYETHYLENE PE
- MONOBLOC STRUCTURE, STAINLESS STEEL AISI-304
- ISO DRAWER FRAMEWORK
- ISO DRAWERS 600 x 400 mm
- ø 100 mm DOUBLE WHEELS, 2 ANTISTATIC, 2 WITH BRAKE

Accessories:

- 26LW(RI): CENTRAL SUPPORT
- 33R: DISPOSABLE CONTAINER SUPPORT
- 50: 13 LITRES BIN

MULTIFUNCTION 300 SERIES

**ISO
640 x 480 mm**

ICU, CARDIOLOGY, DRESSING, OUTPATIENTS CLINICS



STATIC TROLLEYS

DIVIDERS (Optional)



H5/8, (3)H10/4, H20/4



Product code: 3133 W, 365 W.26LW(RI).46.50

Product description:

- DIMENSIONS 640 x 480 mm
- HEIGHT 980 mm (WORK SURFACE)
- WORK SURFACE, POLYETHYLENE PE
- MONOBLOC STRUCTURE, STAINLESS STEEL AISI-304
- ISO DRAWER FRAMEWORK
- ISO DRAWERS 600 x 400 mm
- ø 100 mm DOUBLE WHEELS, 2 ANTISTATIC, 2 WITH BRAKE

Accessories:

- 26LW(RI): CENTRAL SUPPORT
- 46: 3 GLOVE BOX HOLDER
- 50: 13 LITRES BIN

MULTIFUNCTION 300 SERIES



STATIC TROLLEYS

ISO
640 x 480 mm

Product code: 3102 B, 322 B.32R

Product description:

- DIMENSIONS 640 x 480 mm
- HEIGHT 980 mm (WORK SURFACE)
- WORK SURFACE, POLYETHYLENE PE
- MONOBLOC STRUCTURE, STAINLESS STEEL AISI-304
- ISO DRAWER FRAMEWORK
- ISO DRAWERS 600 x 400 mm
- ø 100 mm DOUBLE WHEELS, 2 ANTISTATIC, 2 WITH BRAKE

Accessories:

- 32R: SIDE SUPPORT



Product code: 3107 G, 333 G. 32R

Product description:

- DIMENSIONS 640 x 480 mm
- HEIGHT 980 mm (WORK SURFACE)
- WORK SURFACE, POLYETHYLENE PE
- MONOBLOC STRUCTURE, STAINLESS STEEL AISI-304
- ISO DRAWER FRAMEWORK
- ISO DRAWERS 600 x 400 mm
- á 100 mm DOUBLE WHEELS, 2 ANTISTATIC, 2 WITH BRAKE

Accessories:

- 32R: SIDE SUPPORT



Product code: 3129 W, 365 W.32R

Product description:

- DIMENSIONS 640 x 480 mm
- HEIGHT 980 mm (WORK SURFACE)
- WORK SURFACE, POLYETHYLENE PE
- MONOBLOC STRUCTURE, STAINLESS STEEL AISI-304
- ISO DRAWER FRAMEWORK
- ISO DRAWERS 600 x 400 mm
- ø 100 mm DOUBLE WHEELS, 2 ANTISTATIC, 2 WITH BRAKE

Accessories:

- 32R: SIDE SUPPORT



MULTIFUNCTION 300 SERIES

ISO
900 x 630 mm

TRAUMATOLOGY, DRESSING, SURGERY



Product code: 3850 G, N365 G.26(RA4N.A5).32NR

Product description:

- DIMENSIONS 900 x 630 mm
- HEIGHT 980 mm (WORK SURFACE)
- WORK SURFACE, POLYETHYLENE PE
- MONOBLOC STRUCTURE, STAINLESS STEEL AISI-304
- ISO DRAWER FRAMEWORK
- ISO DRAWERS 400 x 600 mm
- ø 100 mm DOUBLE WHEELS, 2 ANTISTATIC, 2 WITH BRAKE

Accessories:

- 26(RA4N.A5): BINS DISPENSERS
- 32NR: SIDE SUPPORT



Product code: 3820 B, N365 B.26LW(R1N).32NR

Product description:

- DIMENSIONS 900 x 630 mm
- HEIGHT 980 mm (WORK SURFACE)
- WORK SURFACE, POLYETHYLENE PE
- MONOBLOC STRUCTURE, STAINLESS STEEL AISI-304
- ISO DRAWER FRAMEWORK
- ISO DRAWERS 400 x 600 mm
- ø 100 mm DOUBLE WHEELS, 2 ANTISTATIC, 2 WITH BRAKE

Accessories:

- 26LW(R1N): CENTRAL SUPPORT
- 32NR: SIDE SUPPORT

PLASTER TROLLEY 300 SERIES

**ISO
900 x 630 mm**

Stainless steel AISI-304 monobloc structure.

Polyethylene work surface, PE

LXN Stainless steel protection, resists contact with fiberglass bandages



Product code: 3892 W, N322 W.LXN.26(R12N.R2N).37NR.49

Product description:

- DIMENSIONS 900 x 630 mm
- HEIGHT 980 mm (WORK SURFACE)
- WORK SURFACE, POLYETHYLENE PE
- MONOBLOC STRUCTURE, STAINLESS STEEL AISI-304
- ISO DRAWER FRAMEWORK
- ISO DRAWERS 400 x 600 mm
- ø 100 mm DOUBLE WHEELS, 2 ANTISTATIC, 2 WITH BRAKE

Accessories:

- LXN: STAINLESS STEEL PROTECTION
- 26(R12N.R2N): BANDAGE SHELF ROLLS – CENTRAL SHELF
- 37NR: FIXING RAIL
- 49: IPLASTER-BASIN SUPPORT

MRI ROOM

ISO
640 x 480 mm



STATIC TROLLEYS



3 Tesla

Certificate of MRI conformance

DIVIDERS (Optional)



H5/8, (3)H10/4, H20/4



Product code: 3607MRI G, 365MRI G.26(RA4.A5).32R

Product description:

- DIMENSIONS 640 x 480 mm
- HEIGHT 980 mm (WORK SURFACE)
- WORK SURFACE, POLYETHYLENE PE
- MONOBLOC STRUCTURE, STAINLESS STEEL AISI-304
- ISO DRAWER FRAMEWORK
- ISO DRAWERS 600 x 400 mm
- ø 100 mm WHEELS, 4 WITH BRAKE AND ANTISTATIC

Accessories:

- 26(RA4.A5): BINS DISPENSERS
- 32R: SIDE SUPPORT

MRI ROOM



STATIC TROLLEYS

Trolley manufactured with antimagnetic materials
4 antistatic wheels with brake



ISO
640 x 480 mm



3 Tesla

Certificate of MRI conformance

Product code: 3102MRI B, 322MRI B.32R

Product description:

- DIMENSIONS 640 x 480 mm
- HEIGHT 980 mm (WORK SURFACE)
- WORK SURFACE, POLYETHYLENE PE
- MONOBLOC STRUCTURE, STAINLESS STEEL AISI-304
- ISO DRAWER FRAMEWORK
- ISO DRAWERS 600 x 400 mm
- ø 100 mm WHEELS, 4 WITH BRAKE AND ANTISTATIC

Accessories:

- 32R: SIDE SUPPORT

Product code: 3129MRI Y

Product description:

- DIMENSIONS 640 x 480 mm
- HEIGHT 980 mm (WORK SURFACE)
- WORK SURFACE, POLYETHYLENE PE
- MONOBLOC STRUCTURE, STAINLESS STEEL AISI-304
- ISO DRAWER FRAMEWORK
- ISO DRAWERS 600 x 400 mm
- ø 100 mm WHEELS, 4 WITH BRAKE AND ANTISTATIC

Accessories:

- 32R: SIDE SUPPORT

OPERATING ROOM



ANAESTHESIA 300 SERIES

ISO
640 x 480 mm

EMERGENCY ROOM, OPERATING ROOM, REANIMATION AREA



STATIC TROLLEYS



DIVIDERS (Optional)



H5/16, (3)H10/4, H20/4

Product code: 3607 G, 365 G.26(RA4.A5).32R

Product description:

- DIMENSIONS 640 x 480 mm
- HEIGHT 980 mm (WORK SURFACE)
- WORK SURFACE, POLYETHYLENE PE
- MONOBLOC STRUCTURE, STAINLESS STEEL AISI-304
- ISO DRAWER FRAMEWORK
- ISO DRAWERS 600 x 400 mm
- ø 100 mm DOUBLE WHEELS, 2 ANTISTATIC, 2 WITH BRAKE

Accessories:

- 26(RA4.A5): BINS DISPENSERS
- 32R: SIDE SUPPORT

ANAESTHESIA 300 SERIES

OPERATING ROOM

ISO
640 x 480 mm



STATIC TROLLEYS



DIVIDERS (Optional)



H5/16, (3)H10/4, H20/4

Product code: 3684 B, 365 B.26(RA4.A5.R2).30.32L

Product description:

- DIMENSIONS 640 x 480 mm
- HEIGHT 980 mm (WORK SURFACE)
- WORK SURFACE, POLYETHYLENE PE
- MONOBLOC STRUCTURE, STAINLESS STEEL AISI-304
- ISO DRAWER FRAMEWORK
- ISO DRAWERS 600 x 400 mm
- ø 100 mm DOUBLE WHEELS, 2 ANTISTATIC, 2 WITH BRAKE

Accessories:

- 26(RA4.A5.R2): BINS DISPENSERS – CENTRAL SHELF
- 30: FOLD DOWN TRAY
- 32L: SIDE SUPPORT

ANAESTHESIA 300 SERIES

ISO
640 x 480 mm



STATIC TROLLEYS



Product code: 3613 W, 362M6 W.26(RA4.A5).32R

Product description:

- DIMENSIONS 640 x 480 mm
- HEIGHT 980 mm (WORK SURFACE)
- WORK SURFACE, POLYETHYLENE PE
- MONOBLOC STRUCTURE, STAINLESS STEEL AISI-304
- ISO DRAWER FRAMEWORK
- ISO DRAWERS 600 x 400 mm
- ø 100 mm DOUBLE WHEELS, 2 ANTISTATIC, 2 WITH BRAKE

Accessories:

- 26(RA4.A5): BINS DISPENSERS
- 32L: SIDE SUPPORT

EMERGENCY



EMERGENCY SQM SERIES

ISO
480 x 480 mm



EXTENSIVE MOBILITY TROLLEYS



DIVIDERS (Optional)

Product code: Q260 R, SQ161R.12.14.28.32R.34.36.40

Product description:

- DIMENSIONS 480 x 480 mm
- HEIGHT 1015 mm (WORK SURFACE)
- WORK SURFACE, POLYETHYLENE PE
- MONOBLOC STRUCTURE, STAINLESS STEEL AISI-304
- SQUARE DRAWERS
- ELEVATED PUSH HANDLE
- ø 125 mm SIMPLE WHEELS, 2 ANTISTATIC, 2 WITH BRAKE

Accessories:

- 12: MANUAL LOCKING
- 28: ROTATING ARM
- 32R: SIDE SUPPORT
- 34: CARDIAC MASSAGE BOARD
- 36: GAS BOTTLE STAND
- 40: SERUM HOLDER

EMERGENCY SQM SERIES



EXTENSIVE MOBILITY TROLLEYS

Security, manual locking
Numbered disposable seals



SIMPLE WHEELS



Product code: Q272 B, SQ161B.12.14.34.36.40

Product description:

- DIMENSIONS 480 x 480 mm
- HEIGHT 1015 mm (WORK SURFACE)
- WORK SURFACE, POLYETHYLENE PE
- MONOBLOC STRUCTURE, STAINLESS STEEL AISI-304
- SQUARE DRAWERS
- ELEVATED PUSH HANDLE
- ø 125 mm SIMPLE WHEELS, 2 ANTISTATIC, 2 WITH BRAKE

Accessories:

- 12: MANUAL LOCKING
- 34: CARDIAC MASSAGE BOARD
- 36: GAS BOTTLE STAND
- 40: SERUM HOLDER

100 SERIES



EXTENSIVE MOBILITY TROLLEYS

SIMPLE WHEELS ø 125 mm (Optional)

Security, manual locking. Numbered disposable seals



WS12



ISO 640 x 480 mm

Product code: I272 B, B.AXWD12.WS12

Product description:

- DIMENSIONS 640 x 480 mm
- HEIGHT 1015 mm (WORK SURFACE)
- WORK SURFACE, POLYETHYLENE PE
- MONOBLOC STRUCTURE, STAINLESS STEEL AISI-304
- ISO DRAWER FRAMEWORK
- ISO DRAWERS 600 x 400 mm
- ELEVATED PUSH HANDLE
- ø 125 mm SIMPLE WHEELS, 2 ANTISTATIC, 2 WITH BRAKE

Accessories:

- 12: MANUAL LOCKING
- 34: CARDIAC MASSAGE BOARD
- 36: GAS BOTTLE STAND
- 40: SERUM HOLDER

Product code: I208 R, I62M6 R.12.14.28.32R.34.36.40

Product description:

- DIMENSIONS 640 x 480 mm
- HEIGHT 1015 mm (WORK SURFACE)
- WORK SURFACE, POLYETHYLENE PE
- MONOBLOC STRUCTURE, STAINLESS STEEL AISI-304
- ISO DRAWER FRAMEWORK
- ISO DRAWERS 600 x 400 mm
- ELEVATED PUSH HANDLE
- ø 125 mm SIMPLE WHEELS, 2 ANTISTATIC, 2 WITH BRAKE

Accessories:

- 12: MANUAL LOCKING
- 28: ROTATING ARM
- 32R: SIDE SUPPORT
- 34: CARDIAC MASSAGE BOARD
- 36: GAS BOTTLE STAND
- 40: SERUM HOLDER

EMERGENCY 100 SERIES

**ISO
640 x 480 mm**



EXTENSIVE MOBILITY TROLLEYS

DOUBLE WHEELS



DIVIDERS (Optional)



H5/16, (3)H10/4, H20/4



Product code: 1296 B, 165 R, 12, 14, 28, 32R, 34, 36, 40

Product description:

- DIMENSIONS 640 x 480 mm
- HEIGHT 1005 mm (WORK SURFACE)
- WORK SURFACE, POLYETHYLENE PE
- MONOBLOC STRUCTURE, STAINLESS STEEL AISI-304
- ISO DRAWER FRAMEWORK
- ISO DRAWERS 600 x 400 mm
- ELEVATED PUSH HANDLE
- ø 125 mm DOUBLE WHEELS, 2 ANTISTATIC, 2 WITH BRAKE

Accessories:

- 12: MANUAL LOCKING
- 14: NUMBERED DISPOSABLES
- 28: ROTATING ARM
- 32R: SIDE SUPPORT
- 34: CARDIAC MASSAGE BOARD
- 36: GAS BOTTLE STAND
- 40: SERUM HOLDER

TREATMENT



TREATMENT SQM SERIES

**ISO
480 x 480 mm**

HOSPITALISATION UNITS



EXTENSIVE MOBILITY TROLLEYS

SIMPLE WHEELS



Product code: Q102 W, SQ122 W.32R.50

Product description:

- DIMENSIONS 480 x 480 mm
- HEIGHT 1015 mm (WORK SURFACE)
- WORK SURFACE, POLYETHYLENE PE
- MONOBLOC STRUCTURE, STAINLESS STEEL AISI-304
- SQUARE DRAWERS
- ELEVATED PUSH HANDLE
- ø 125 mm SIMPLE WHEELS, 2 ANTISTATIC, 2 WITH BRAKE

Accessories:

- 32R: SIDE32R: SIDE SUPPORT
- 50: 13 LITRES BIN SUPPORT
- 50: 13 LITRES BIN



Product code: Q106 W, SQ161 W.32R.46.50

Product description:

- DIMENSIONS 480 x 480 mm
- HEIGHT 1015 mm (WORK SURFACE)
- WORK SURFACE, POLYETHYLENE PE
- MONOBLOC STRUCTURE, STAINLESS STEEL AISI-304
- SQUARE DRAWERS
- ELEVATED PUSH HANDLE
- ø 125 mm SIMPLE WHEELS, 2 ANTISTATIC, 2 WITH BRAKE

Accessories:

- 32R: SIDE SUPPORT
- 46: 3 GLOVES BOX HOLDER
- 50: 13 LITRES BIN



Product code: Q136 B, SQ161 B.26LW(RIQ).33R.46.50

Product description:

- DIMENSIONS 480 x 480 mm
- HEIGHT 1015 mm (WORK SURFACE)
- WORK SURFACE, POLYETHYLENE PE
- MONOBLOC STRUCTURE, STAINLESS STEEL AISI-304
- SQUARE DRAWERS
- ELEVATED PUSH HANDLE
- ø 125 mm SIMPLE WHEELS, 2 ANTISTATIC, 2 WITH BRAKE

Accessories:

- 26LW(RIQ): CENTRAL SUPPORT
- 33R: DISPOSABLE CONTAINER SUPPORT
- 46: 3 GLOVES BOX HOLDER
- 50: 13 LITRES BIN

TREATMENT 100 SERIES

ISO
640 x 480 mm

HOSPITALIZATION UNITS



EXTENSIVE MOBILITY TROLLEYS

DOUBLE WHEELS



DIVIDERS (Optional)



H5/8, (3)H10/4, H20/4

Product code: I606 B, I65 B.32R.46.50

Product description:

- DIMENSIONS 640 x 480 mm
- HEIGHT 1005 mm (WORK SURFACE)
- WORK SURFACE, POLYETHYLENE PE
- MONOBLOC STRUCTURE, STAINLESS STEEL AISI-304
- ISO DRAWER FRAMEWORK
- ISO DRAWERS 600 x 400 mm
- ELEVATED PUSH HANDLE
- ø 125 mm DOUBLE WHEELS, 2 ANTISTATIC, 2 WITH BRAKE

Accessories:

- 32R: SIDE SUPPORT
- 46: 3 GLOVES BOX HOLDER
- 50: 13 LITRES BIN50: 13 LITRES BIN

UNIT DOSE MEDICATION 100 SERIES

**ISO
640 x 480 mm**



EXTENSIVE MOBILITY TROLLEYS



DIVIDERS (Optional)
SEPARADORES (Opcional)



H10/4, H20/4

Product code: I809W, I12M6 W.32R.46.50

Product description:

- DIMENSIONS 640 x 480 mm
- HEIGHT 1005 mm (WORK SURFACE)
- WORK SURFACE, POLYETHYLENE PE
- MONOBLOC STRUCTURE, STAINLESS STEEL AISI-304
- ISO DRAWER FRAMEWORK
- ISO DRAWERS 600 x 400 mm
- CASSETTE FOR MEDICATION, INTERCHANGEABLE, 5 PATIENTS
- KEY LOCKING
- ELEVATED PUSH HANDLE
- ø 125 mm DOUBLE WHEELS, 2 ANTISTATIC, 2 WITH BRAKE

Accessories:

- 32R: SIDE SUPPORT
- 46: 3 GLOVES BOX HOLDER
- 50: 13 LITRES BIN

UNIT DOSE MEDICATION 100 SERIES

ISO
640 x 480 mm




Product code: I894 B, I14M6 B.32R.46.50

Product description:

- DIMENSIONS 640 x 480 mm
- HEIGHT 1005 mm (WORK SURFACE)
- WORK SURFACE, POLYETHYLENE PE
- MONOBLOC STRUCTURE, STAINLESS STEEL AISI-304
- ISO DRAWER FRAMEWORK
- ISO DRAWERS 600 x 400 mm
- CASSETTE FOR MEDICATION, INTERCHANGEABLE, 5 PATIENTS
- KEY LOCKING
- ELEVATED PUSH HANDLE
- ø 125 mm DOUBLE WHEELS, 2 ANTISTATIC, 2 WITH BRAKE

Accessories:

- 32R: SIDE SUPPORT
- 46: 3 GLOVES BOX HOLDER
- 50: 13 LITRES BIN



3 Seymour Court, Tudor Road, Manor Park, Runcorn, Cheshire WA7 1SY
01928 571 801 | info@distinctivemedical.com | www.distinctivemedical.com

Registered in the UK. Registered office: 3 Seymour Court, Tudor Road, Manor Park, Runcorn, Cheshire WA7 1SY
VAT Registration No. 775518887