



distinctive
medical

Tamper Evident Collection



ORDERING FROM DISTINCTIVE MEDICAL

Ordering from Distinctive Medical
could not be easier....

Website www.distinctivemedical.com

Phone 01928 571801

Email info@distinctivemedical.com

Mail Distinctive Medical Ltd
3 Seymour Court
Tudor Road
Manor Park
Runcorn
Cheshire, WA7 1SY

Payment

To ensure prompt action for your order or enquiry,
please quote your purchase order number, name and
phone number on all orders and correspondence.

We accept BACS, CHAPS and corporate or personal
debit and credit cards. Please do not send cash through
the post.



Please note: We do not accept American Express or
Electron as forms of payment.

Pricing

All our prices exclude VAT and delivery. We are
constantly striving to improve our product quality and
reserve the right to amend any products and pricing
shown in this catalogue.

ISO 9001

Distinctive Medical meets the standards and requirements
set by ISO 9001:2015



Delivery options

We aim to maintain competitive everyday prices,
we do not have a minimum order value. We do however,
have a delivery charge of £9.95 per order.

All delivery charges are exclusive of VAT.

Standard parcel

Subject to stock availability we aim to deliver your order
the next working day following receipt of your order*
All orders are subject to a flat shipping rate of £9.95.

Premium parcel services

Next day before 12:00 noon

Subject to a delivery charge of £30.00 per order**

Next day before 10:00am

Subject to a delivery charge of £35.00 per order**

Palletised or Bulky Items

Palletised or bulky items are shipped on a 2-3 day
service and are charged at £49.50 per pallet.
Next day services are available on such deliveries
and there is a £18.00 per pallet surcharge for this service.

*Orders must be placed by 15:30 for delivery the
following day.

**For premium parcel services orders must be
placed by 14:00 to meet timed deliveries.

The above prices and times are applicable to
Mainland UK customers. Overseas customers,
including Channel Islands and Irish Republic, should
contact Customer Services for carriage charges.

Please note: Working day excludes Saturday, Sunday
and Bank Holidays.

Returns Policy

You may return any product to us within 7 days of
receipt of delivery providing the product is in a
resaleable condition. If an item is faulty or incorrect
then we must be notified in writing within 24 hours of
receipt of delivery. Customers are responsible for the
cost of returning a non-faulty product to us. We
reserve the right to refuse credit or replacement
goods if goods are returned damaged or of a
sub-standard condition. Customised products and
SPL product codes are non-returnable.



Tamper Evident Bags

TAMPER EVIDENT RANGE

Pharmacy Gusset

Product SKU: GTPHARM

Tamper resistant bag for the transportation of items from the pharmacy. Lasts for over 2000 cycles of use. The Envopak Pharmacy Gusset is designed to safely transport and store drugs when in transit from the pharmacy to the hospital ward.

When used with tamper evident tags this becomes a tamper resistant bag and indicates to users if there has been unauthorised access. The bag features a yellow modular chamber where the zip on the bag slots in. The zip can then be secured by either a PS2 or a PS4 Envopolyseal. 1. (See page 7 & 8)



Pharmacy Carrier

Product SKU: CRPHARM

Tamper-evident lockable bag for pharmacy use. 330 denier PVC coated polyester. Can be reused up to 2000 times. The pharmacy carrier is a lockable bag that comes in green with the message 'return to the pharmacy department' printed on the front. This bag has been given a durability rating of heavy duty and can be used an estimated 2000 times.

This pharmacy carrier also comes with handles for ease of use. The inside of the bag is padded and has inserts for the safe storage of drugs. This secures the items from unnecessary breakages and spillages of its contents. Making the bag tamper resistant is simple, all you have to do is pull the zip across to the yellow modular chamber and slot the zip in. Then add a PS2 or PS4 Envopolyseal to secure.



Medical Record Carrier

Product SKU: MCMED

Medical record security tamper evident bag. 550 denier PVC coated polyester. 2000 cycles of reusability. This medical record security tamper evident bag deters unauthorised access with a sealed lock. The bag can be secured using either our PS2 or PS4 Envopolyseals – working in adherence to the Caldicott Report.

If the seal is broken it indicates that someone who isn't permitted has opened the bag and its contents could be tampered with. This bag is designed to store documents such as medical records. With handles on



TAMPER EVIDENT RANGE

Cytotoxic Pharmacy Carrier

Product SKU: CRCYTOX

Tamper evident when used with seals. 332 denier PVC coated polyester. Can be used for up to 2000 cycles. The Cytotoxic Carrier tamper bag safely transports cytotoxic drugs to the hospital ward from the pharmacy.

The bag is yellow making it easily recognisable that its use is for cytotoxic drugs only. This is in line with the ISOPP Standards of Practice Transport of Cytotoxic guidelines. Features the message 'cytotoxic handle with care return to the pharmacy' printed on it in black.



A3 Flat Tamper Evident Pouch

Product SKU: A3Flat

Tamper evident security bag for documents. Photo-sized window to display information. Can be reused up to 2000 times.

Available in blue or red this modular design Envopak has been designed to be a security bag for documents.

With a secure closure for the zip that can be locked in place with a tamper-evident seal. This allows you to be confident when sending important documents.

This bag works with PS2 and PS4 seals that come numbered.



A4 Flat Tamper Evident Pouch

Product SKU: A4Flat

Tamper evident locking security bag. Photo-sized window for personalisation. Lasts up to 2000 cycles. This A4 locking security bag gives its contents added security whilst in transit.

This is due to the tamper-evident seals that lock the zip in place. The seals that work with this bag are the PS2 and PS4 polyseals.

They are numbered to provide an additional check to ensure contents haven't been tampered with when travelling between locations.



TAMPER EVIDENT RANGE

A5 Flat Tamper Evident Pouch

Product SKU:A5Flat

Tamper evident bag security. Photo-sized window on the outside to display information. Lasts 2000 cycles. Safely transport documents from A to B with the A5 flat. The light-duty PVC provides a protective barrier to the contents of the bag whilst remaining light and easy to store when not in use.

The bag features a yellow modular chamber where the zip can be locked into place with a tamper-evident seal. The seals that work with this bag are the PS2 and PS4 polyseals. The tamper-evident bag security allows the recipient of the contents to be made aware if there has been any unauthorised access. The seals that work with this bag are numbered for additional security.



A6 Flat Tamper Evident Pouch

Product SKU:A6Flat

Tamper evident lockable pouch. Information window. Lasts for up to 2000 uses. The A6 modular flat is a lockable pouch that provides a strong barrier for the storage of documents and other items whilst they are in transit. It is available in three colour options – blue, red and clear. Please tell us your choice when ordering.

This bag can be made tamper evident and is compatible with our PS2 and PS4 polyseals. These seals are numbered so ideal if you need to provide an audit trail for the bag's contents.



A7 Flat Tamper Evident Pouch

Product SKU:A7Flat

Tamper evident security pouch. Plastic window to display information. Lasts up to 2000 uses. This A7 modular flat is designed to contain small-sized items such as cash and receipts. Due to its dimensions, this bag also works with airflow tube systems.

Ensure the contents have not been tampered with whilst they have been in transit with the use of the PS2 and PS4 polyseals that are compatible with this bag.



TAMPER EVIDENT RANGE

Key Wallet Flat

Product SKU: IKWI04A

Tamper evident lock bag. 550 denier PVC coated polyester. Can be reused up to 2000 times. Available in blue and plastolene – for visibility of contents. This lock bag is designed to be a secure storage solution for keys.

On the outside of the bag, there is a window where information can be displayed. When the zip is locked with a PS2 or PS4 seal in the yellow chamber it is made tamper evident. The seals are numbered for additional checking and security.



Fridge Item Enclosed Return To Pharmacy Bag

Product SKU: /

See our Fridge Item Enclosed Return To Pharmacy Bag. These bags read “Fridge Item Enclosed. Return To Pharmacy Production Unit”. The bag has a modular window, a yellow PS4 seal and an orange nylon zip slot.

Bag Specifications:
Size: 483mm x 381mm
Material: Light Duty
Colour: Orange




PS2 Polyseals

Product SKU: PS2 + Colour

The PS2 Envo polyseals are used to secure the Envopak bags, to protect its contents during storage and transportation. The security seals are packed in boxes of 1,000 and available in Blue, Green, Red, Orange, White and Yellow – Please select preferred colour when ordering.

The security seals are printed with 7 digit sequential numbers for additional security. The number being reordered and tracked allows for the construction of an audit trail.



TAMPER EVIDENT RANGE

PS4 Polyseal

Product SKU: PS4

The PS4 Envopolyseals are used to secure the Envopak bags, to protect its contents during storage and transportation. The PS4 has multiple benefits compared to the PS2 security seal: It has a newer design for easier application and removal from the new yellow chamber

Same price as the PS2 but it is easier to use. The numbers and size of the seal is bigger, making it easier to read which may improve the security level of the seal.


The PS4 security seals are packed in boxes of 1,000 and available in Blue, Green, Red, Orange, White and Yellow – Please select your preferred colour when ordering. The security seals are printed with 7-digit sequential numbers for additional security. The number being reordered and tracked allows for the construction of an audit trail.



AlertVoid Security Bags

Product SKU(s): I0440 | I0442 | I0444 | I0472 | I0476 | I0473
| I0475 | I0477

These Security Bags are available in numerous sizes and in packs of 25 or rolls of 500. They are also available in clear, grey, or white. The dual serial numbers provide traceability and chain-of-custody reporting to control and track prescriptions, samples, and other items needing security during transport. The Bags are closed with a permanent pressure-sensitive adhesive and contain a large write-on area.



Keep Safe Bags

Product SKU: / KTSE 1, KTSE 2, KTSE 3

These tamper-proof plastic bags feature identification codes and barcodes. Each bag has a tear-off receipt with a code and barcode, and a VOID message appears if the bag has been tampered with. Tamper-proof plastic bags provide security against a range of attacks including moisture, solvent and mechanical.

Available in white or clear plastic each bag features its own unique identification code and barcode. This means its contents can have an audit trail whilst they are in transit. These bags are single-use and feature a tear-off receipt with a code and barcode on for an accurate reference. Available in quantities of 500.



TAMPER EVIDENT RANGE

High Alert Tamper Indicating Bag

Product SKU: 17501

Available 100 per package, these bags are great for error prevention initiatives! Clear bags with unique red trim and a white write-on area.

Tamper evident feature protects the integrity of high-alert medications.

The bags have a perforated tear-off receipt tag for documentation and are self-sealing.



Self-Sealing Tamper-Evident Bags

Product SKU(s): 17812 | 17499 | 1714 | 1718 | 1720 | 7524 | 7526

Tamper-evident feature protects the integrity of medications. Perforated tear-off receipt tag for documentation. Contains a white write-on area.

Thickness – 1.5-mil

Closure – Self-Sealing

Sizes Available – 381 x 76mm (17812) 95 x 159mm (7514) 95 x 343mm (7518) 64 x 267mm (17499) 152 x 254mm (7524) 165 x 167mm (7526)



Self-Sealing Tamper-Evident Bags (Amber)

Product SKU: 17859

These amber bags provide you with UV light protection. The tamper-evident bags contain a white write-on area. Perforated tear-off receipt.

Thickness – 1.5 mil. Closure – Self-Sealing

Quantity – 100 per pack

Dimensions – 95 x 159mm



TAMPER EVIDENT RANGE

Compact Cold Chain Vaccine Transportation Bag 3 Litre Vaccine Bag

Product SKU: 2990

The Compact Cold Chain Vaccine Transportation Bag is a medicine and vaccine carrier that is ideal for transporting cold chain items short distances such as between wards or clinics. This insulated bag is designed for the transportation of chilled products and extending the time between 2-8°C. With a new design, the vaccine bag designs bring improvements on durability, infection control properties, insulation and handling.

Internal Dimensions: W355 x D245 x H60mm

External Dimensions: W370 x D270 x H110mm

Insulation Thickness: 45mm

Estimated capacity: 20 HPV packs – without the use of gel packs/coolants



10 litre Small Cold Chain /Vaccine Transportation Bag

Product SKU: 2991

A lightweight medical cool bag designed for the transportation of chilled products and extending the time between 2-8°C. With a new design, the vaccine bag designs bring improvements on durability, infection control properties, insulation, and handling.

Internal Dimensions: W395 x D255 x H170mm

External Dimensions: W410 x D280 x H220mm

Insulation Thickness: 45mm

Estimated Capacity: 80 HPV bags – without the use of gel packs or coolants.



20 litre Large Cold Chain /Vaccine Transportation Bag

Product SKU: 2992 / 299T - WITH TROLLEY OPTION

The Large Cold Chain/Vaccine Transportation Bag, is designed for transporting chilled products and extending the time between 2-8°C. These medium-sized 20-litre sized bags are ideal for transporting cold chain items for longer periods of time, such as between wards or clinics or for short community visits.

- Internal Dimensions: W395 x D255 x H290mm
- External Dimensions: W410 x D280 x H340mm
- Insulation Thickness: 45mm
- Estimated Capacity: 120 HPV packs – without the use of gel packs and coolants.



TAMPER EVIDENT RANGE

30 Litre Extra Large Cold Chain /Vaccine Transportation Bag

Product SKU: 2993 / 2993T - with trolley

Our 30 Litre Extra Large Cold Chain /Vaccine Transportation Bag is our biggest yet and is ideal for large quantities. This carrier is ideal for transporting cold chain items for a longer period of time, such as for day-long community visits. These larger bags can be heavy when fully loaded with coolant and chilled contents, thus an integrated trolley option is also available.

- Internal Dimensions: W415 x D315 x H300mm
- External Dimensions: W430 x D340 x H350mm
- Insulation Thickness: 45mm
- Estimated Capacity: 220HPV packs – without use of gel packs or coolants
- Trolley Size: 360 x 300mm (please specify if you would like this option when ordering)



Medication Transport Case

Product SKU: 19964

- Compact – Soft-sided medication carrier protects temperature-sensitive contents during transport. When used as directed, the case is proven to maintain 2-8°C for up to 12 hours based on testing in the ambient temperature range of 10°-21°C.
- Reliable Protection – Includes four PharmaChill™ cold packs and cold pack layering instructions for maximum drug stability. Can be secured with security seals in zipper pulls.
- Dimensions – 265W x 216H x 216mmD



Transport Cooler

Product SKU: 20185 /20186

Our Transport Cooler is available in 2 sizes. Good In, Bad Out – Confidently transport temperature-sensitive meds without fear of patient safety compliance failures or leaks. Sturdy latches and a leak-resistant gasket provide a tight seal that keeps cold in and prevents spills from escaping out. Durable – Transport often results in wear and tear that can compromise your container, but this heavy-duty plastic shell is tough to handle.

Dimensions –

Small (product code 20185) – 6.5l – 338 x 282 x 220mm

Large (product code 20186) – 15l – 420 x 325 x 284mm



TAMPER EVIDENT RANGE

Lined Transport Bag

Product SKU: 19731 / 19732 / 19733

Our Lined Transport Bag is available in 3 sizes.

- Safer Transport – Insulated transport tote features a tamper-evident lock that allows users to safely transport temperature-sensitive medication. Tamper-evident seals (not included) can be purchased separately.
- Document Protection – The front pouch allows for storage of important documents.
- Easy to Clean – Liner attaches with Velcro®, removes easily for cleaning.
- Secure – with PS2 or PS4 tamper evident seals

Dimensions –

Small (product code 19731) – 273 x 203 x 209mm

Medium (product code 19732) – 276 x 178 x 308mm

Large (product code 19733) – 330 x 222 x 308mm



Insulated Transport Bag

– Blood, Specimens and Human Tissue – UN3373 approved.

Product SKU: EB04.003

COOL'S is an isothermal bag that allows the transport in an organized way and under optimal conditions of diagnostic samples, from the extraction centre to the laboratory. The bag is certified according to "UN 3373", that enables the transport of biological substances class B. Optionally EB04.003 Cool's bag can be delivered with 6 or 12 pieces of reusable cooling gel packs EB09.007 (25,5 x 10,5 cm)

Size: 42x26.5x9.5cm. –3 compartments for test tube racks. Each one has a capacity for 2 racks (100 test tubes), with a size of 26.5x14x27 cm. All detachable compartments have a handle, identifier and sponge to absorb spills.

- Measures: 44 x 29 x 39 cm
- Capacity: 49,76 L
- Weight: 5.30 kg
- Material: Tarpaulin. Water resistant. Washable
- Colour: Blue





Tamper Evident Seals

TAMPER EVIDENT RANGE

Pull Fly Baby Seal – Consecutively Numbered

Product SKU: 7924 / 7928 / 7917 / 7923 / 7926 / 7925

Fits smaller boxes and areas with limited space where other types of seals may have difficulty sealing or breaking for removal. Average breaking strength – 10-1/2 lb. Available on NHS supply chain.

Material – Polypropylene plastic

Quantity – 100 per package

Tools needed for removal – None.

Autoclavable – Up to 325°F (163°C)

Also available in Blue, Orange, Red, White and Yellow.




Easy Pull-Tight Seals – Numbered

Product SKU: 18319

One-piece, pull-tight, the adjustable seal is designed to enhance safety by showing evidence of tampering once the seal has been removed. The self-locking fastener is designed to be easily applied and removed.

Requires no tools for removal, saving valuable time when accessing medication or supplies. All-plastic security tag features an engineered break point that allows quick access in emergencies. The thin sealing strap fits most locking rings and equipment latches. Less waste! Produced in strips of 10 with no plastic waste when separated from the strip. Available in: Black, Blue, Green, Red, White, Yellow



Easy Pull-Tight Seals – Unnumbered

Product SKU: 18320

One-piece, pull-tight, the adjustable seal is designed to enhance safety by showing evidence of tampering once the seal has been removed. The self-locking fastener is designed to be easily applied and removed. Requires no tools for removal, saving valuable time when accessing medication or supplies.

All-plastic security tag features an engineered break point that allows quick access in emergencies. The thin sealing strap fits most locking rings and equipment latches. Less waste! Produced in strips of 10 with no plastic waste when separated from the strip. Available in : Black, Blue, Green, Red, White, Yellow



TAMPER EVIDENT RANGE

Envotite Un-numbered Seals

Product SKU: PESB / PESY / PESR / PESG

The Envopak pull tight unnumbered seal provides tamper-evident protection for trolleys, boxes, cases, fire extinguishers and first aid boxes. The security seal contains 1,000 seals per box and is available in Yellow, Red, Green, and Blue – Please select your preferred colour when ordering.

This economical and simple seal has a low breaking strength which has been designed to be easily removed by hand.



Envotite Numbered Seals

Product SKU: NSYDM | NSRDM | NSBDM | NDWDM | NSGDM

The Envotite Numbered Seals provide tamper-evident protection for trolleys, boxes, cases, fire extinguishers and first aid boxes. The security seal contains 1,000 seals per box. Available in Blue, Green, Red, White, and Yellow – Please select your preferred colour when ordering. Low breaking strength, easy to break by hand.

The tamper-evident security seal has large clear consecutive numbers which facilitates maintenance checking and identification of tampering or unauthorised access.



Handy Padlock Seal – Numbered

Product SKU: I8317B / I8317G / I8317K / I8317R / I8317W / I8317Y

These seals can be used to secure tote boxes and trolleys. Remove by hand. Allows quick and easy access. Available consecutively numbered and unnumbered.

Autoclavable up to 446°F (230°C)

100 seals per package.

Available in Black (I8317K), Blue (I8317B), Green (I8317G), Red (I8317R), White (I8317W) or Yellow (I8317Y).



TAMPER EVIDENT RANGE

Handy Padlock Seal – Un-Numbered

Product SKU: 18318B / 18318G/ 18318K / 18318R / 18318W / 18318Y

These seals can be used to secure tote boxes and trolleys. Remove by hand. Allows quick and easy access. Available consecutively numbered and unnumbered.

Autoclavable up to 446°F (230°C)

100 seals per package.

Available in Black (18318K), Blue (18318B), Green (18318G), Red (18318R), White (18318W) or Yellow (18318Y).




Spring Lok Pre-Printed Locking Tag

Product SKU: 7906 / 7912 / 7914 / 7901 / 7910 / 7903

Label and secure medical supplies in seconds with the Spring Lok Pre-Printed Locking Tag a tamper evident security seal. These seals are ideal for many medical departments with the following messages printed onto the seal: Return to OR, Floor Tray, Broken, Discard, Dirty, Clean, Loaner, Quick Turnover, Pull Set for Service, Repair and Sharpen.

These messages allow for the quick identification of what items need maintenance.

No additional tools are needed to break the lock. The tamper-evident security seal is made from polypropylene and is autoclavable – up to 163°C. Store each seal with ease using the slit doughnut stringing hole.



MiniFlex Seal

Product SKU: MiniFlex

All plastic pull-up seals are suitable for securing medical trolleys and carts, cases, curtain buckles and tote boxes. A smooth pull-up design allows for speedy application. The seals have large clear numbers for easy identification.

The tear line allows for removal by hand. Thin strap ideal for small apertures/ eyelets. Contains no metal parts – is recyclable.



TAMPER EVIDENT RANGE

Ringlok Seal

Product SKU: R215G / R215Y (AVAILABLE IN GREEN (R215G) AND YELLOW (R215Y))

These seals can be used for sealing fire doors & vehicle doors. Strong fixed-length seal. The seals are easy to apply and remove by hand. Alpha characters on the strap to prevent cutting and re-gluing. Laser marked for maximum security.



Grip Seal

Product SKU: GRIPSEAL

The Gripseal security seals are ideal for securing clinical waste or high-security items. This security seal is available in boxes of 1,000 in red, yellow and blue – please select your preferred colour when ordering.

Easy grip tear line and easy feed label holder. Tail keeper to hold each tooth with no reverse slipping.

Chamber holds multiple teeth to eliminate “wriggle room”. Centralised locking chamber design for more efficient packing.



Omega II Seal

Product SKU:

The Omega II security seal has a simple pull-up design which allows for speedy application. This security seal is available in boxes of 1,000 in white and yellow – please select your preferred colour when ordering. Removable by hand as the seal has a tear line.

The tamper-evident security seal has large clear numbers which facilitate maintenance checking and identification of tampering or unauthorised access.



TAMPER EVIDENT RANGE

Pull Tight Unnumbered Seal (pack of 100)

Product SKU: 7827 / 7825 / 7931 / 7821 / 7828 / 7823

Excellent seal for code carts and emergency kits that need security with quick entry. Average breaking strength – 11-14 lb. Printable area – One side of the seal tab is textured to permit the writing of verification information (e.g., date, time, initials).

Material – Polypropylene plastic

Quantity – 100 per package

Tools needed for removal – None

Autoclavable – Up to 325°F (163°C)

Available in blue (7827), green (7825), purple (7931), red (7821), white (7828), yellow (7823)



Pull Tight Numbered Seal (pack of 100)

Product SKU: 7816 / 7817 / 7930 / 7814 / 7813 / 7930

Excellent seal for code carts and emergency kits that need security with quick entry. Average breaking strength – 11-14 lb. Printable area – One side of the seal tab is textured to permit the writing of verification information (e.g., date, time, initials).

Material – Polypropylene plastic

Quantity – 100 per package

Tools needed for removal – None

Autoclavable – Up to 325°F (163°C)



Pull-Tight Padlock Seal

Product SKU:

These seals combine the benefits of both the “Pull Tight” and “Padlock” types of seals. Average breaking strength – 12-1/2 lb. Printable area – Seals can be written on with a permanent marker to identify contents, date, time, initials, etc. Unnumbered.

Material – Polypropylene plastic

Quantity – 100 per package

Tools needed for removal – None

Autoclavable – Up to 243°F (117°C)



TAMPER EVIDENT RANGE

Secur-Pull Breakable Seal – Consecutively Numbered or Numbered

Product SKU: 7909-10 / 7908-10 / 8188-10 / 7907-10 / 8187-10

Consecutively numbered seals. Available in black (7909-10), blue (7908-10), green (8188-10), red (7907-10) and yellow (8187-10). Package of 100. These tamper-proof security seals have an average breaking strength of 24lb.

A breaking point at the base of the pull tab ensures no tools are needed to break the seal but can be done by hand. The seals are consecutively numbered allowing for audit trails to be created for the items they are securing. The tamper-proof security seals are autoclavable up to 163°C and are made using polypropylene plastic.



Ring-Pull II Bi-Directional Seal

Product SKU: 8141 / 8140 / 8143 / 8144 / 8145 / 8146 / 8147

Consecutively numbered breakable security seals. These seals are available in black (8141), blue (8140), green (8143), orange (8144), red (8145), white (8146) and yellow (8147). They are available in packages of 100 seals. The insert on these breakable security seals is bi-directional meaning the seal can be secured from either side. This can speed up the process of securing items as no need to reposition the seal to an enter-only side.

The seal can be broken by hand with the rods attached to the tail tearing away. The flat surface on the seal can be written on using a permanent marker. This is useful if the date, time or initials need to be recorded.



Ring-Pull Seal

Product SKU: 8127 / 8123 / 8122 / 8124 / 8120 / 8125 / 8121

Consecutively numbered seals. Available in blue (8123), black (8127), yellow (8121), green (8122), orange (8124), red (8120) and white (8125). Packages of 100 seals. These tamper-indicating seals can be quickly broken due to their ring pull lock. This makes the seals ideal for emergency equipment where speedy access is vital. With a double wall feature that protects the seal when locked into place on the inside. The outer wall protects tampering from the outside.

The average breaking strength is 23lb meaning no tools are needed to break the seal. Made from polypropylene plastic the tamper indicating seals are autoclavable up to 163°C. The seals come consecutively numbered in a choice of six colours making them ideal for a range of organisational systems.



TAMPER EVIDENT RANGE

Secur-Grip Security Seal – Numbered

Product SKU: 7851-10 / 7842-10 / 8190-10 / 7841-10 / 8189-10

Numbered seals. Available in black (7851-10), blue (7842-10), green (8190-10), red (7841-10) and yellow (8189-100). Package of 100 seals. These heavy-duty plastic seals have an average breaking strength of 120 lbs. **This means scissors or a knife is needed for removal.** Made from a strong polypropylene plastic with pronged tails. The strength of these seals and printed identification numbers makes them an effective tool against tampering.

Length of seal tail: 273mm

Total length of seal: 343mm

Information such as date, time and initials can be written on the surface of the seal in permanent marker. These heavy-duty plastic seals are autoclavable up to 163°C.



Extra-Large Heavy-Duty Padlock Seals – Numbered

Product SKU: 8331 / 8330 / 8332 / 8329 / 8333 / 8334

Average breaking strength – 20 lbs. These seals are made from Acetal.

Tools are needed for the removal of these seals. They are autoclavable – Up to 171°C

Colours – Black (8331), Blue (8330), Green (8332), Red (8329), White(8333) and Yellow(8334)

Quantity – 100 per package

Dimensions – 35mm x 52mm



Foil Tamper Indicating Vial Seal

Product SKU: 7866, 7867, 7887

These foil vial seals offer flexibility when sealing a vial. This ensures there is a safe and strong moisture, particulate and bacterial barrier when appropriately applied.

The foil vial seals come with an adhesive on them to allow for application before freezing or refrigeration. Available in a choice of colours to suit any colour-coded system such as days of the week and type of medication. Choose from blue, green, orange, red and yellow. Please let us know what colour you would like when ordering.

Available in quantities of 1000 per roll in a dispensing box.





Tamper Evident Tape and labels

TAMPER EVIDENT RANGE

Tamper Tape Dispenser & Tamper Resistant Tape

Product SKU: 8004 | 8000 | 8006 | 8002

Easily store and dispense up to two rolls of tamper-evident tape in a single device for quick, convenient access anywhere in your facility. Eliminates aggravating searches for loose tape rolls when securing specimen vials, medical containers and other critical care items that need to be safely sealed.


Tape available to purchase separately in Red SKU:8000 | Clear SKU:8002



Permanent Tamper-Evident Labels

Product SKU: 7957B / 7957R / 7958B / 7958R

Our range of permanent tamper-evident labels is designed for any number of applications including the securing of primary and secondary packaging. Available in a range of sizes and colours. The permanent tamper-evident labels leave an "OPEN VOID" message on the surface, once removed. If you require a tamper-evident label that does not leave any message on the surface, our non-residue tamper-evident labels are the ideal solution.



Non-Residue Tamper Evident Labels

Product SKU: 5995, 5997, 5999

Our range of non-residue, tamper-evident labels are designed to be used on hospital furniture and equipment where you do not want to leave a void message on the item itself. Instead, the void message is immediately displayed on the label once removed. If you are looking for tamper-evident labels that leave a message on the surface, our permanent tamper-evident labels are the ideal solution. Available in a range of colours and sizes. Our tamper-evident labels can also be customised to your specification by adding personalised text or logo on the labels.

Pack size 1000

5995	85 x 25	Red or Blue
5997	55 x 20	Red or Blue
5999	140 x 38	Red or Blue



TAMPER EVIDENT RANGE

Tamper Evident Tape Strips

Product SKU: 12281

Use with boxes, specimen vials and other hard-to-secure items to help deter and show tampering.

Tape moulds easily to almost any item. Red tape with Black print. Tape reads: Integrity Seal, Sealed by and Date.

Dimensions: 81mm x 35mm
Available 100 strips per box.



Syringe Seal

Product SKU: 8001

- For use on all syringes – oral and injectable.
- Seals can also be used on vials, boxes, otic/ophthalmic containers, etc.
- Paper seals mould easily to the syringe cap or sheath.
- Wrinkle and tear when removed.
- Non-sterile
- Seal reads: DO NOT USE IF SEAL IS BROKEN
- Dimensions: 89mm x 13mm
- 1,000 seals per roll in a dispensing box.



Tamper Indicating Seal, Preprinted

Product SKU: 18813

Tamper Indicating Seals can be used to cover and protect all types of syringes, oral dispensers, boxes, vials, jars and more. Pre-printed. Yellow warning seal with Black print reads "Do not use if seal is broken".

These seals are also available without text. (SKU:18812). The Tamper Indicating Seals deter and show signs of tampering by tearing when removed from drug packaging.



TAMPER EVIDENT RANGE

Sterile Chemo Bag Port Seal

Product SKU 18695

Safe, sterile, tamper-evident seals for admixture ports on IV bags, bottles, prefilled syringes and vials. Protect against contamination and dosing errors.

Alerts clinicians when a drug has been added or modified. Single-use, removal leaves a printed 'OPENED' warning.



Red Bag Port Seal

Product SKU: 18296R

Sterile, tamper-evident Foil bag port seals protect the point of entry of IV bags containing chemotherapy drugs. Pharmacy safety seals with pre-printed message cover the IV bag port to prevent medication contamination and accidental double dosing.

The tamper-indicating foil outer layer provides visible evidence of a needle puncture and makes it easy to see the bag is free from tampering. Removal of the seal's outer layer leaves a printed "OPENED" warning and cannot be reapplied. Seals maintain adhesion when refrigerated or frozen in cold drug storage conditions.
1,000 seals per roll.



Blue Bag Port Seal

Product SKU: 18296B

Sterile, tamper-evident Foil bag port seals protect the point of entry of IV bags containing chemotherapy drugs. Pharmacy safety seals with pre-printed message cover the IV bag port to prevent medication contamination and accidental double dosing.

The tamper-indicating foil outer layer provides visible evidence of a needle puncture and makes it easy to see the bag is free from tampering. Removal of the seal's outer layer leaves a printed "OPENED" warning and cannot be reapplied. Seals maintain adhesion when refrigerated or frozen in cold drug storage conditions.



TAMPER EVIDENT RANGE

Belly Button Bag Port Seal

Product SKU: ST6101

Safe, sterile Green bag port seal covers the belly button port of an IV bag, providing tamper indication and protection from contamination and accidental double dosing.

The foil outer layer provides visible evidence of a needle puncture to see the bag is free from tampering, ensuring the admixture's integrity from preparation to administration. (1,000 seals per roll)



Self-Adhesive Test Tube Seal

Product SKU: 7963B / 7963R

Wrap over test tube caps or stoppers for visual evidence of tampering to lab samples and specimens. Specifically designed for test tube use. Maintain tube integrity seals act as a deterrent to tampering.

Tear-off receipt can be applied to records and packaging for documentation and help ensure a correct match amongst groups of tubes.

Consecutive numbering and consecutive number barcodes provide traceability for samples in tubes. Scan barcodes to simplify tracking and eliminate handwritten control logs. Fits Tubes from 20-25mm.

Self-Adhesive Test Tube Seal Dimensions: 140mm x 19mm
250 per dispensing box.





Tamper Resistant Caps

TAMPER EVIDENT RANGE

Friendly and Safe Vials with Snap Caps Or Child Resistant Caps - 9 Dram – 60 Dram

Product SKU: Various codes

Order vials and caps separately or with the caps already attached to each vial. Quantity - 100 vials with caps per package.

Caps and Vials can be made tamper evident with tamper resistant tape.



Tamper Resistant Cap

Product SKU: I0294

Available 250 per package. The tamper-resistant cap has a ring that breaks away from the cap when it is opened. One size fits all of the Amber plastic bottles. (Bottles sold separately) The Amber plastic bottles that accompany these caps are made of polyethylene terephthalate (PET) plastic.



Tamper Resistant Cap, Black

Product SKU: I0068

This Tamper Resistant Cap, Black is available 100 per case. The tamper-resistant cap comes in black it has an aluminium foil seal inside.

The ring breaks off to show it has been opened. The Amber plastic bottles that accompany these caps are made of polyethylene terephthalate (PET) plastic.



Tamper Resistant Cap, White

Product SKU: I0067

This Tamper Resistant Cap, Black is available 100 per case. The tamper-resistant cap comes in white it has an aluminium foil seal inside. The ring breaks off to show it has been opened. The Amber plastic bottles that accompany these caps are made of polyethylene terephthalate (PET) plastic.



TAMPER EVIDENT RANGE

Tamper Evident Cap

Product SKU: 7716

Tamper-evident safety caps, made of three-piece aluminium construction, deter tampering. Once the cap is crimped on the vial, the contents cannot be accessed without removing the cap and leaving evidence of access.

Vials sold separately (#7707, 7709, 7721, 7726 and 7731). Our unit dose packaging inventory includes a variety of oral dispensers, vials and caps to satisfy all of your liquid unit dose packaging needs from start to finish. 5,000 caps per package.



Child Resistant Screw Cap for Cylinder Plastic Bottles

Product SKU: 10412

Our Child Resistant Screw Cap for Cylinder Plastic Bottles are available 100 per pack. Using the Child Resistant Screw Cap helps to provide safety closure for plastic bottles. One size of the Child Resistant Screw Cap perfectly fits all bottle sizes. Bottles available on our website.



Sterile Tamper-Evident Caps for Luer Lock Syringes

Product SKU: 17361R

Sterile Tamper-Evident Caps for Luer Lock Syringes. Better Integrity – Sterile Tamper-Evident Tip Cap helps maintain the integrity of a luer lock syringe by restricting access to the tip and leaving behind a plastic sleeve when the cap is removed to indicate use or tampering. Convenient – Individually sterile packaged for small quantity preparations. Compatible – Caps fit all luer lock syringes.

Colour – Red


Quantity – 100 per package – sku: 17361

25 per package – sku: 31-20

Colour – White

Quantity – 100 per package – sku: 18955

25 per package – sku: 31-50



TAMPER EVIDENT RANGE

Product	Code	Page
Double Glove Dispenser – Stainless Steel	7473-01	28
Easy Empty Crusher Bag®	12679	6
Easy-Grip Tablet Crusher	7336	49
Enema Bottles	Various codes	30
Ergo-Mat Anti-Fatigue Mat	5289-01	42
Face Mask Dispenser	Various codes	28
Foam Sealing Tray 7006/Budget/Compact	6468	39
Foam Sealing Tray for #7024-7026 Deep Round-Cavity Blisters	6464	39
Foam Sealing Tray for 6 Down X-Large	6471	40
Foam Sealing Tray For 7 Down	6469	37
Foam Sealing Tray for Extra High Barrier Blister	6472	37
Foam Sealing Tray for Extra High-Barrier Blister Block	6460	38
Foam Sealing Tray for High Barrier S-M-L	6470	38
Foam Tray For Shallow Round Blisters	6461	38
Foil Label for 7021-7026 Round-Cavity Extra High-Barrier Blister	6125	39
Foil Label for High Barrier X-Large	6122	38
Foil Label High Barrier Blister Laser	6118	39
Friendly and Safe Vial (9 dram – 60 dram)	Various codes	42
Gel Packs – 1 litre	GELPACKS	9
Glass Bottles	Various codes	30
Glass Ointment Slab, 1/4 Inch Thick	3098	15
Glass Ointment Slab, 3/4 Inch Thick	3413	15
Glass Stirring Rod, 10 Inch	3018	15
Glass Stirring Rod, 12 Inch	3009	16
Glass Stirring Rod, 8 Inch	3017	16
Glassine Paper	3016	16
Glassine Sheets	3097	16
Glove Box Dispenser Double – Space Saver	GS-108	29
Hinged Door Locking Refrigerator Box, Small, Single Lock	3746-01	11
Horizontal Refrigerator/Freezer Thermometer	16088	13
Insulated Transport Bag	EB04.003	7
Insulated Transport Box – Large	19807	9
Insulated Transport Box – Small	19808	9
IV Tubing cover	SPL-26461	4
IV Cover – Light Amber	7590L	4

TAMPER EVIDENT RANGE

Product	Code	Page
IV Cover with Slit at Sealed End – Dark Amber	19644	4
IV Cover with Slit at Sealed End – Green	7579	5
IV Cover with Slit at Sealed End – Light Amber	Various codes	5
IV Tubing Cover, Amber, 72 In. Long	7617	4
Jumbo 7 Day AM/PM Med Planner	1132C	21
Kango® High Polyurethane Seat Chair without Tilt with Footrest	8423K	43
Kango® Low Polyurethane Seat Chair With Tilt and Casters	8422K	42
Large Weekly Medication Organizer	7128-01	21
Laser Foil Labels for Blister Pack	6126	36
Laser Foil Labels for Blisters	6121	37
Latex Free Tourniquets On A Roll	10311	45
Lid for Narrow Graduated Medication Cups, 400 Pkg.	5160	32
Lined Transport Bag	Various codes	9
Long Probe Thermometer	17844	13
Long Tamper Syringe Seal – Do Not Use If Seal Is Broken	7998	45
Medium Amber Extra High-Barrier Blister	7014	40
Medication Crushing Cup Set	7337-01	49
Medication Transport Case	19964	8
Memory Monitoring Refrigerator/Freezer Thermometer	10368	12
Monthly Medication Organiser	7109	21
Mortar Card	3035-50	15
Narrow Graduated Medication Cup, Blue, 400 Pkg.	5162	32
Pharmacy Carrier	CRPHARM	26
Paper Label for Compact Blisters Pinfeed	7030	36
Paper Souffle Cup, 1 oz.	7369	31
Paper Souffle Cup, 1/2 oz.	7368	31
Personal Tablet Crusher With Cap, Blue	7335	49
Plastic Bottles	Various codes	30
Plastic Tablet Crusher	7305B	48
Plastic Utility Box, X-Large, 14x2x9	1488	20
Porcelain Mortar and Pestle Set (3 Sizes)	Various codes	49
Portable IV Pole, Large	1575-01	18
Press In Adaptor	Various codes	41
PS2 Polyseals	PS2 + Colour	25
PS4 Polyseal	PS4 + Colour	26
Punch Card Divider	6215-01	36

About Us.

Let's face it, healthcare demands are tougher than ever. Budgets are tight, regulations are even tighter, and patients are counting on you to make the right choices for their care. It's a heavy burden to bear, but **Distinctive Medical** is ready to help.


For more than 30 years, we have been an industry leader specialising in manufacturing, packaging, and distributing unique and hard-to-find healthcare products. We specialise in small package quantities with no minimum order requirements, fast shipping, and hassle-free returns, because we always strive to meet the needs of our customers!

Our website is packed full of innovative new products and practical solutions to everyday problems.

We also offer lots of great ideas for making your routine easier, your results more effective and your job more enjoyable. And if there's a product you need but can't find, contact our Resource Department to find it for you.

Our friendly, knowledgeable Customer Service staff are available from 9:00 a.m. to 5:00 p.m., Monday to Friday. So, give us a look and give us a call on 01928 571801.

You are our number one priority!



Discover the full range of Distinctive Medical products...



**Avalo Medical Trolley
Series Brochure**



**Medical Labels
Brochure**



**The Pharmacy
Collection Brochure**



**Sepsis Storage
Solutions Brochure**